

Test Health Exchange Dashboard (HED) Portal

November 2009
Document version 1.0

© 2009 Compuware Corporation
All Rights Reserved
Confidential and Proprietary

Table Of Contents

Before Using Test Health Exchange Dashboard (HED) Portal.....	1
About the Test HED Portal	3
Workflow for Using Test HED Portal	5
List Test Case Tab.....	7
Add/Edit Test Case Tab.....	9
Test Log Tab	11
Create New Test Case.....	13
Initiate a Test Case	15
Delete a Test Case Response	17
Download a Response Document.....	19
Edit a Response Document	21
How to Edit Test Case Values.....	23

Before Using Test Health Exchange Dashboard (HED) Portal

The following items are required before using Test HED Portal:

- you have prerequisite knowledge of messaging; understanding sent and received payload information
- you have been granted access to Test HED Portal
- you have the appropriate datasets in order to create/edit test cases
- you have identified the custom artifact request and response format of the targeted SP to be simulated by Test HED Portal
- Covisint has completed the backend setup with your Service Providers (SPs)
- a Test HED Portal instance has been established

About the Test HED Portal

Test HED Portal (Test HED) is a custom application used to validate connections between Service Providers (SPs) and the associated Health Exchange Dashboard (HED) portlets on the portal. It consists of an administrative portlet which enables users to define key combinations to create test cases, and to submit test artifact datasets which simulate PHI artifacts coming from certain (SPs) which are then displayed in the appropriate HED portlet.

The integrated multi-application setup includes:

- Test HED Portal instance
- Record Locator Service (RLS)
- Master Patient Index Service (MPI),
- Covisint Connection Administration (CCA)
- Covisint Connect Messaging Hub (Connect).

Test Case Name	Practice Name	Patient ID	Patient Last Name	Provider Name	Artifact Type	Bulk Case	Delete Case
BLAS104H	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104B	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104C	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104D	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104E	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104F	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104G	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104H	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104I	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104J	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104K	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104L	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104M	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104N	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104O	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104P	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104Q	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104R	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104S	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104T	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104U	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104V	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104W	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104X	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104Y	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		
BLAS104Z	4132374-HealthPartners	41000000000000000000	ANDERSON	4142380-Blue Cross Blue Sh...	SLIDELITY		

View larger image >>>

Workflow for Using Test HED Portal

A tester should have a prior knowledge of the custom artifact request and response format of the targeted SP to be simulated by Test HED Portal .

- create a valid response document with the expected results (outside of Test HED Portal system)
- **create a test case**, upload the response document, and provide a set of keys to make it possible to match the response with the artifact request. These keys include:
 - o practice ID (usually a TP_ID number for the practice generated in Connect while patient sync setup is performed),
 - o patient ID (something known to the MPI),
 - o destination service provider ID (pre-configured in Connect)
 - o artifact type (pre-configured in Test HED Portal)
 - o It is possible to use two additional keys for the eligibility checks and claim statuses in the absence of the member patient ID - first and last name of the patient.
- open the HED search portlet

The selected artifact portlet will call RLS. The RLS will:

- create the canonical artifact request document
- lookup the predefined list of SPs for the simulated practice or healthcare community
- distribute the artifact request to the listed simulated SPs by calling Connect synchronous processing engine.

The Covisint Connect application will:

- forward the canonical artifact request to the special synchronous application channel that simulates the SPs.
- Test HED application channel responder parses the canonical artifact request
- retrieves all required keys (see above)
- tries to match them against the existing database records with the configured test cases.

If the test case is found, the preconfigured artifact response will be returned back to Connect engine and later to RLS.

If the test case is not found in the Test HED database, the error response will be created for RLS and error entry will show up in the Test HED responder log.

The RLS will:

- consolidate all of the received responses into a one composite response

Test Health Exchange Dashboard (HED) Portal

- return it back to the appropriate HED portlet.

The portlet will display the patient results on the HED screen accordingly.

List Test Case Tab

List Test Cases								Add/Edit Test Case	Test Log
Practice Name: ALL		Patient ID: ALL	Patient Last Name: ALL	Provider Name: ALL	Artifact Type: ALL				
Page 1 of 16		1 - 10 / 159							
Test Case Name	Practice Name	Patient ID	Patient Last Name	Provider Name	Artifact Type	Edit Case	Delete Case		
ELUGTQH	4132374-HealthPartners	KHO699156031		4142390-Blue Cross Blue SH...	ELIOLITY				
ELUGJA9	4132374-HealthPartners	D1002202103	ANDERSON	4122397-UCare	ELIOLITY				
ELUGJA8	4132374-HealthPartners	D1002202103	ANDERSON	4142390-Blue Cross Blue SH...	ELIOLITY				
ELUGAS7	4132374-HealthPartners	D1006202104	SWANSON	4122397-UCare	ELIOLITY				
ELUGM6	4132374-HealthPartners	D1002201101	JOHNSON	4142390-Blue Cross Blue SH...	ELIOLITY				
ELUGBOGB	4132374-HealthPartners	D100600057042	GRANGER	4122397-UCare	ELIOLITY				
- DEMO_MEDHST5	76111-State of Minnesota De...	205101	Granger	76111-State of Minnesota De...	MEDHST				
- DEMO_MEDHST3	4122395-Medica	207101	Granger	4122395-Medica	MEDHST				
TCBG05	4132374-HealthPartners	D100511157042	GRANGER	4122395-Medica	ELIOLITY				
TCBG04	4132374-HealthPartners	D100411157042	GRANGER	4132374-HealthPartners	ELIOLITY				
Page 1 of 16		1 - 10 / 159							

The List Test Cases tab displays all currently available test cases showing the case key combinations and management action icons. This tab consists of:

- **Sorting** function, to allow for sorting of the test cases by the:
 - source practice ID
 - patient ID
 - SP ID
 - artifact type
- **Edit Case** action icon - opens the case details on the "Add/Edit Test Case" screen, where the user can modify the ID values and/or upload a new response payload for the test case.
- **Delete Case** action icon - deletes the test case from the system. (Delete function is permanent. There is no recycle bin for deleted cases).
- **Filter Options** - The filter drop down options contain only the practice name, patient id, provider name and artifact type.

Add/Edit Test Case Tab

The Add/Edit Test Case tab serves two purposes:

1. Add new test cases by clicking on the tab.
2. Edit existing test cases by clicking the **Edit** icon of the test case from the **List Test Case** tab. When clicked, the edit screen displays the Add/Edit tab with fields populated by data from the selected test case.

Add Test Case

Add Test Case Screen

Test Case Name:
Patient ID:
Patient First Name:
Patient Last Name:
Response Name:
Response Body: [Browse...](#)
Portlet Instance ID: QA
Practice Name/ID:
Provider Name/ID:
Artifact Type:
[Save](#) [Cancel](#)

Required fields on this tab include:

- Test Case Name,
- Either the Patient Id, or the combined fields of Last Name + First Name
- Practice Name
- Provider Name
- Artifact Type

Edit Test Case

Edit Test Case Screen

Test Case Name: ELIGTOH
Patient ID: KH0899156031
Patient First Name:
Patient Last Name:
Response Name:
Response Body: [Browse...](#)
Portlet Instance ID: QA
Practice Name/ID: 41323744HealthPartners
Provider Name/ID: 4142380-Blue Cross Blue Shield of Minn
Artifact Type: ELIGIBILITY
[Save](#) [Cancel](#)

Manage Test Case Responses:

Active	Response Name	Upload Response	Download Response	Delete Response
1	ELIGTOH			

[Add New...](#)
[Save](#) [Cancel](#)


In the Manage Test Case Responses section of the Edit view:

- **Modify** ID values
- **View payloads** are displayed, and only one can be one identified as active at any time

Test Health Exchange Dashboard (HED) Portal

- **Upload response** option - enables the user to select a new test case to replace the existing
- **Download response** option - enables user to view the message response
- **Delete response** option - enables user to delete the message response from the system
- **Add New...** option - enables the user to upload a new test case


Test Log Tab

List Test Cases Add/Edit Test Case Test Log		
Viewing Test Log Since: <input type="text"/>		
Page: 1 of 2149 		
Timestamp	Message ID	Log Entry
10/09/2009 13:10:53.000	8042000230	SUCCESS: Case "M4-BO_DEMO_SRC3 (M4-BO_SOURCES.dmf)" returned for keys SRC:4132374-PRV:4132374-ART...
10/09/2009 13:10:53.000	8042000232	SUCCESS: Case "BO_MEDHST2 (EMPTY)" returned for keys SRC:4122395-PRV:4122395-ART:MEHST-PATD...
10/09/2009 13:10:53.000	8042000231	SUCCESS: Case "BO_MEDHST4 (EMPTY)" returned for keys SRC:4122395-PRV:4122395-ART:MEHST-PATD...
10/09/2009 13:10:53.000	8042000234	SUCCESS: Case "BO_CCD3 (CCD-BO_FULL1)" returned for keys SRC:4382389-PRV:4382389-ART:CCD-PATD:10...
10/09/2009 13:10:53.000	8042000233	SUCCESS: Case "M4-BO_DEMO_SRC5 (M4-BO_SOURCES.dmf)" returned for keys SRC:4142380-PRV:4142380-ART...
10/09/2009 13:10:53.000	8042000229	SUCCESS: Case "BO_MEDHST5 (EMPTY)" returned for keys SRC:76111-PRV:76111-ART:MEHST-PATD:100...
10/09/2009 13:10:53.000	8042000228	SUCCESS: Case "BO_CCD1 (BO_CCD1)" returned for keys SRC:4382389-PRV:4382389-ART:CCD-PATD:1000000...
10/09/2009 13:10:51.000	8042000226	SUCCESS: Case "BO_MMJNT (BO_MMJNT)" returned for keys SRC:4382387-PRV:4382387-ART:MMJNT-PATD:100...
10/09/2009 13:09:59.000	8042000223	SUCCESS: Case "BO_LAB1 (LAB-SHORT)" returned for keys SRC:4372386-PRV:4372386-ART:LABS-PATD:4372...
10/09/2009 13:09:59.000	8042000222	SUCCESS: Case "BO_LAB2 (LAB_BO_BASIC_PANEL.dmf)" returned for keys SRC:4382386-PRV:4382386-ART:LA...

The Test Log tab displays Test HED responder log, used for troubleshooting purposes. Here the user is able to view:

- a particular log in Test HED by precise timestamp
- a particular log in TEST HED Connect message ID
- whether a particular test was successful or failed
- information necessary to debug cases with the wrong combination of keys defined

Using the Test Log Tab:

- Click  to refresh the log with the latest tests

A log entry contains all information needed to set up a test case; practice id, provider id, patient id, artifact type and patient name

A log will show all requests for artifacts based on the sources the patient is associated with in the MPI.

If the log shows SUCCESS, and yet the portlet does not display the data, then there is likely a problem with the data (or less likely a problem with the rls/connect setup). In this case, the Connect msgid can be taken from the log and forwarded to Covisint Messaging Team for investigation.

An error does not necessarily indicate a problem, but rather a test case for that key combination has not been set up.

If the test harness log shows ERROR, it simply means that the patient is associated with a practice/provider for which there is no test case set up.

Create New Test Case

1. Click **Add/Edit Test Case** tab.

The screenshot shows the 'Add/Edit Test Case' form. It includes tabs for 'List Test Cases', 'Add/Edit Test Case', and 'Test Log'. The 'Add/Edit Test Case' tab is selected. The form contains the following fields and controls:

- Test Case Name: Text input field
- Patient ID: Text input field
- Patient First Name: Text input field
- Patient Last Name: Text input field
- Response Name: Text input field
- Response Body: Text input field with a 'Browse...' button
- Portlet Instance ID: Dropdown menu (selected: QA)
- Practice Name/ID: Dropdown menu
- Provider Name/ID: Dropdown menu
- Artifact Type: Dropdown menu
- Save: Button
- Cancel: Button

2. Key in all fields, as all are required.
3. In the *Response Name* field, create a name for the response that is easily recognizable by you.
4. In the *Response Body* field, click **Browse...** and select the response document you wish to attach to this case. (Response documents are created outside of Test HED Portal).
5. From the drop down menu, select:
 - a. Practice Name / ID
 - b. Provider Name / ID
 - c. Artifact Type
6. Click **Save**. The screen refreshes and the test case is created.

RESULT:

You have successfully created a test case.

Initiate a Test Case

Complete the following steps in order to validate a request and response in Test HED Portal .

Look up a patient (via administrative or clinical views tab of your portal, depending on the artifact type)

View the appropriate artifact portlet - portlet will state 'x of x responded'.

The first number here is the number of simulated sources responded to the artifact request (i.e., how many test cases matched) and the second number is the number of sources providing this artifact type that the patient is associated with.

Delete a Test Case Response

1. Navigate to the **List Test Case** tab.

List Test Cases Add/Edit Test Case Test Log							
Practice Name: ALL Patient ID: ALL Patient Last Name: ALL Provider Name: ALL Artifact Type: ALL							
Page 1 of 16							
Test Case Name	Practice Name	Patient ID	Patient Last Name	Provider Name	Artifact Type	Edit Case	Delete Case
ELJGTOH	4132374-HealthPartners	410099156031		4142380-Blue Cross Blue SH...	ELIGIBILITY		
ELJGJAB	4132374-HealthPartners	41002202103	ANDERSON	4122387-UCare	ELIGIBILITY		
ELJGJAB	4132374-HealthPartners	41002202103	ANDERSON	4142380-Blue Cross Blue SH...	ELIGIBILITY		
ELJGAS7	4132374-HealthPartners	41006202104	SWANSON	4122387-UCare	ELIGIBILITY		
ELJGMA6	4132374-HealthPartners	41002201101	JOHNSON	4142380-Blue Cross Blue SH...	ELIGIBILITY		
ELJGBO02B	4132374-HealthPartners	4100600057042	GRANGER	4122387-UCare	ELIGIBILITY		
- DEMO_MEDHST5	76111-State of Minnesota De...	205101	Granger	76111-State of Minnesota De...	MEDHST		
- DEMO_MEDHST3	4122395-Medica	207101	Granger	4122395-Medica	MEDHST		
TCB006	4132374-HealthPartners	4100511157042	GRANGER	4122395-Medica	ELIGIBILITY		
TCB004	4132374-HealthPartners	4100411157042	GRANGER	4132374-HealthPartners	ELIGIBILITY		

2. Click (edit icon) in the Edit column of the test case you wish to edit. The test case is displayed in the Add/Edit Test Case screen with fields populated by data from the selected test case.

List Test Cases Add/Edit Test Case Test Log			
Add/Edit Test Case			
Test Case Name: ELJGTOH	Portfolio Instance ID: QA		
Patient ID: 410099156031	Practice Name/ID: 4132374-HealthPartners		
Patient First Name:	Provider Name/ID: 4142380-Blue Cross Blue Shield of Minn...		
Patient Last Name:	Artifact Type: ELIGIBILITY		
Manage Test Case Responses:			
Active	Response Name	Upload Response	Delete Response
1	ELJGTOH		
Save Cancel			

3. Click in the Delete Response column for the selected response name. The Select file pop up box is displayed.
4. Click **Yes** to confirm deletion.
5. Click **Save**.

RESULT:

You have successfully deleted a test case response.

Download a Response Document

1. Navigate to the **List Test Case** tab.

List Test Cases

Add/Edit Test Case

Test Log

Practitioner Name: ALL

Patient ID: ALL

Patient Last Name: ALL

Provider Name: ALL

Artifact Type: ALL


Page 1 of 16

1 - 10 / 159

Test Case Name	Practice Name	Patient ID	Patient Last Name	Provider Name	Artifact Type	Edit Case	Delete Case
ELUG00A	4123274-HealthPartners	010009956031	ANDERSON	4123280-Blue Cross Blue SH...	ELUGILITY		
ELUGA9	4123274-HealthPartners	010000220103	ANDERSON	4122287-UCare	ELUGILITY		
ELUGA8	4123274-HealthPartners	010000220104	ANDERSON	4122380-Blue Cross Blue SH...	ELUGILITY		
ELUGA7	4123274-HealthPartners	010006220104	SHAWSON	4122287-UCare	ELUGILITY		
ELUGA6	4123274-HealthPartners	010000220101	JOHNSON	4122380-Blue Cross Blue SH...	ELUGILITY		
ELUG00B	4123274-HealthPartners	01000000507042	GRANDER	4122287-UCare	ELUGILITY		
_EMMO_MEDHST5	76111-State of Minnesota De...	209101	Orange	76111-State of Minnesota De...	MEDHST		
_EMMO_MEDHST3	4122395-Medica	070101	Orange	4122395-Medica	MEDHST		
TC8006	4123274-HealthPartners	0100011157042	GRANDER	4122395-Medica	ELUGILITY		
TC8004	4123274-HealthPartners	01000411157042	GRANDER	4123274-HealthPartners	ELUGILITY		

Page 1 of 16

1 - 10 / 159

2. Click  (edit icon) in the Edit column of the test case you wish to edit. The test case is displayed in the Add/Edit Test Case screen with fields populated by data from the selected test case.

[List Test Cases](#)
[Add/Edit Test Case](#)
[Test Log](#)

[Add/Edit Test Case](#)

Test Case Name:

Patient ID:

Patient First Name:

Patient Last Name:

Portlet Instance ID:

Practice Name/ID:

Provider Name/ID:

Artefact Type:


Manage Test Case Responses:

Add New...

Active	Response Name	Upload Response	Download Response	Delete Response
1 <input type="radio"/>	ELJGTOGH			

Save

Cancel

- Click  in the Download Response column for the selected response name.
- Click **Open**. The response is displayed in a new window.

RESULT:

You have successfully downloaded a response document.

Edit a Response Document

1. Navigate to the **List Test Case** tab.

List Test Cases								Add/Edit Test Case		Test Log		
Practice Name: ALL		Patient ID: ALL		Patient Last Name: ALL		Provider Name: ALL		Artefact Type: ALL				
Page 1 of 16												1 - 10 / 159
Test Case Name	Practice Name	Patient ID	Patient Last Name	Provider Name	Artefact Type	Edit Case	Delete Case					
ELJGTOH	4132374-HealthPartners	KH0899156031		4142380-Blue Cross Blue SH...	ELJGILITY							
ELJGJAB	4132374-HealthPartners	E1002202103	ANDERSON	4122387-UCare	ELJGILITY							
ELJGJAB	4132374-HealthPartners	E1002202103	ANDERSON	4142380-Blue Cross Blue SH...	ELJGILITY							
ELJGAS7	4132374-HealthPartners	E1006202104	SWANSON	4122387-UCare	ELJGILITY							
ELJGMA6	4132374-HealthPartners	E1002201101	JOHNSON	4142380-Blue Cross Blue SH...	ELJGILITY							
ELJGBO02B	4132374-HealthPartners	E100600057042	GRANGER	4122387-UCare	ELJGILITY							
- DEMO_MEDHST5	78111-State of Minnesota De...	205101	Granger	78111-State of Minnesota De...	MEDHST							
- DEMO_MEDHST3	4122395-Medica	207101	Granger	4122395-Medica	MEDHST							
TCB006	4132374-HealthPartners	E100511157042	GRANGER	4122395-Medica	ELJGILITY							
TCB004	4132374-HealthPartners	E100411157042	GRANGER	4132374-HealthPartners	ELJGILITY							
Page 1 of 16												1 - 10 / 159

2. Click (edit icon) in the Edit column of the test case you wish to edit. The test case is displayed in the Add/Edit Test Case screen with fields populated by data from the selected test case.

List Test Cases	Add/Edit Test Case	Test Log										
<p>Add/Edit Test Case</p> <p>Test Case Name: ELJGTOH</p> <p>Patient ID: KH0899156031</p> <p>Patient First Name:</p> <p>Patient Last Name:</p> <p>Portlet Instance ID: QA</p> <p>Practice Name/ID: 4132374-HealthPartners</p> <p>Provider Name/ID: 4142380-Blue Cross Blue Shield of Minn...</p> <p>Artefact Type: ELJGILITY</p> <p>Add Name...</p>												
<p>Manage Test Case Responses:</p> <table border="1"> <thead> <tr> <th>Active</th> <th>Response Name</th> <th>Upload Response</th> <th>Download Response</th> <th>Delete Response</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>ELJGTOH</td> <td></td> <td></td> <td></td> </tr> </tbody> </table> <p>Save Cancel</p>			Active	Response Name	Upload Response	Download Response	Delete Response	1	ELJGTOH			
Active	Response Name	Upload Response	Download Response	Delete Response								
1	ELJGTOH											

3. Click in the Upload Response column for the selected response name. The Select file pop up box is displayed.
4. Click **Browse...** and select the edited response document you wish to use as replacement.
5. Click **Save**. The pop up box is closed.
6. Click **Save**.

RESULT:

You have successfully edited a response document.

How to Edit Test Case Values

1. Navigate to the **List Test Case** tab.

Test Case Name	Practice Name	Patient ID	Patient Last Name	Provider Name	Artifact Type	Edit Case	Delete Case
ELJGTOH	4132374-HealthPartners	4132374-HealthPartners	ANDERSON	4142380-Blue Cross Blue SH...	ELJGBLTY		
ELJGJA9	4132374-HealthPartners	4132374-HealthPartners	ANDERSON	4142380-Blue Cross Blue SH...	ELJGBLTY		
ELJGJA8	4132374-HealthPartners	4132374-HealthPartners	ANDERSON	4142380-Blue Cross Blue SH...	ELJGBLTY		
ELJGAS7	4132374-HealthPartners	4132374-HealthPartners	SWANSON	4142380-Blue Cross Blue SH...	ELJGBLTY		
ELJGMA6	4132374-HealthPartners	4132374-HealthPartners	JOHNSON	4142380-Blue Cross Blue SH...	ELJGBLTY		
ELJGBO02B	4132374-HealthPartners	4132374-HealthPartners	GRANGER	4142380-Blue Cross Blue SH...	ELJGBLTY		
-DEMO_MEDHST5	76111-State of Minnesota De...	205101	Granger	76111-State of Minnesota De...	MEDHST		
-DEMO_MEDHST3	4122395-Medica	207101	Granger	4122395-Medica	MEDHST		
TCB006	4132374-HealthPartners	4100511157042	GRANGER	4122395-Medica	ELJGBLTY		
TCB004	4132374-HealthPartners	4100411157042	GRANGER	4132374-HealthPartners	ELJGBLTY		

2. Click (edit icon) in the Edit column of the test case you wish to edit. The test case is displayed in the Add/Edit Test Case screen with fields populated by data from the selected test case.

Add/Edit Test Case

Test Case Name: ELJGTOH
 Patient ID: 4132374-HealthPartners
 Patient First Name:
 Patient Last Name:
 Portlet Instance ID: QA
 Practice Name/ID: 4132374-HealthPartners
 Provider Name/ID: 4142380-Blue Cross Blue Shield of Minn...
 Artifact Type: ELJGBLTY

Manage Test Case Responses:

Active	Response Name	Upload Response	Download Response	Delete Response
1	ELJGTOH			

Save Cancel

3. Perform one or more of the following:

IF YOU WISH TO MODIFY THE...	THEN:
Patient id number	key in the correct number in the Patient ID field.
Test case name	key in the correct name in the Test Case Name field.
Patient name	key in the correct patient name in the first name and/or last name field.
Portlet instance ID	key in the correct Instance ID in the Portlet Instance ID field.
Practice Name / ID	select a new Practice Name / ID from the drop down menu.
Provider Name / ID	select a new Provider Name / ID from the drop down menu.
Artifact type	select a new artifact type from the drop down menu.

4. Click **Save**.